

Wrist Worn Fall Detector



Introduction

The wrist worn fall detector provides reassurance when it's needed the most. It can be used as a wrist worn pendant in addition to its fall detection capabilities.

Using intelligent algorithms, the device constantly monitors the activities of the individual and automatically raises an alert should it believe that a fall has taken place.

The fall detector has been designed for user comfort with contoured ergonomics and elegant aesthetics.



Features and Benefits

- Intelligent fall detection algorithm.
- Low profile design.
- Water resistant.
- User replaceable battery.
- Pendant call function.
- Traditional watch strap fastening.
- Red and Green reassurance LEDs.
- Bi-directional TAO technology.
- Digital Heartbeat.
- Daily battery self-test and low battery reporting.

Compatibility

Compatible with TAO compatible alarm receivers

Reach IP	Carer Response	Advent XT2 Digital Grouped Living Solutions
----------	----------------	---

Legrand Assisted Living and Healthcare

Blyth Workspace, Commissioners Quay, Quay Road, Blyth, Northumberland NE24 3AF.

Tel: 01670 352 371 | Fax: 01670 362 807 | Email: sales@tynetec.co.uk | Web: www.tynetec.co.uk

Specification

Power	• 1 year battery life (Typical) • User replaceable CR2032 Lithium Battery
Conformity	• CE Marked
Standards	Compliant with the following standards: • Electromagnetic Compatibility (EMC) Emissions & Immunity Standards: • EN 61000-6-3:2007+A1:2011 • EN 50130-4:2011+A11:2014 • Radio Equipment Directive: ETSI EN 300 220-3-1 V2.1.1 (2016-12) • Social Alarm System Standard: EN 50134-3:2012
Telecare Transceiver	• 869.2125MHz Class 1.5
IP Rating	• Water Resistant
Operating Temperature	• Typical 15°C to +35°C
Warranty	• 2 years

Dimensions

Dimensions	• 49mm (h) x 30mm (w) x 9mm (d)
Weight	• 29 grams

Product Codes

- **ZXT805** - Wrist Worn Fall Detector

As our policy is one of constant product improvement the right is therefore reserved to modify product specifications without prior notice.

Legrand Assisted Living and Healthcare

Blyth Workspace, Commissioners Quay, Quay Road, Blyth, Northumberland NE24 3AF.

Tel: 01670 352 371 | Fax: 01670 362 807 | Email: sales@tynetec.co.uk | Web: www.tynetec.co.uk