


Door Contact



Introduction

Designed to monitor entry and exit activity in both individual homes and grouped living environments

When a monitored door is opened, it transmits an alarm and the LED will illuminate red for 4 seconds. When the alarm is acknowledged by the Telecare system the LED will illuminate green for 2 seconds. A new alarm will not be triggered until the door has been closed and opened again.



Features and Benefits

- Fixed Trigger device.
- Used to support the care of dementia patients.
- Integrated left and right-side reed switches.
- Surface magnetic door contact.
- Red alarm transmit LED indicator.
- Green alarm received LED Indicator.
- In/Out passage delay option.
- Daily battery self test and low battery reporting.
- Digital Heartbeat.
- Bi-directional TAO technology.

Compatibility

Compatible with TAO compatible alarm receivers

Reach IP	Carer Response	Advent XT2 Digital Grouped Living Solutions
----------	----------------	---

Legrand Care

Innovation House, Silvertown Court, Northumberland Business Park, Cramlington, NE23 7RY.

Tel: 01670 352 371 | Email: sales@tynetec.co.uk | Web: www.tynetec.co.uk

Registered Office: Great King Street North, Birmingham, B19 2LF | Registered in England No. 00115834

Version No:1.2 Door Contact. 09/23

Specification

Power	• 2 No. AAA Alkaline batteries
Conformity	• CE Marked
Standards	Compliant with the following standards: • Electromagnetic Compatibility (EMC) Emissions & Immunity Standards: • EN 61000-6-3:2007+A1:2011 • EN 50130-4:2011+A11:2014 • Radio Equipment Directive: ETSI EN 300 220-3-1 V2.1.1 (2016-12) • Social Alarm System Standard: EN 50134-3:2012
Telecare Transceiver	• 869.2125MHz Class 1.5
Operating Temperature	• +5°C to +40°C
Warranty	• 2 years

Dimensions

Dimensions	• 85mm (h) x 34mm (w) x 24mm (d)
Weight	• 55 grams

Product Codes

- **ZXT830** - Door Contact

As our policy is one of constant product improvement the right is therefore reserved to modify product specifications without prior notice.

Legrand Care

Innovation House, Silverton Court, Northumberland Business Park, Cramlington, NE23 7RY.

Tel: **01670 352 371** | Email: sales@tynetec.co.uk | Web: www.tynetec.co.uk

Registered Office: Great King Street North, Birmingham, B19 2LF | Registered in England No. 00115834

Version No:1.2 Door Contact. 09/23